


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	02/08/20	Interface Bias	: +4.5kV
Time Processed :	07:27:45 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	GP (33-41)	B.conc	: 50%H ₂ O/50%MeOH
Cat No. :	RP20091		
Lot No. :	E0319920B		
Theoretical MW :	1016.18		
Observed MW :	1016.0		