


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	12/14/18	Interface Bias	: +4.5kV
Time Processed :	11:27:38 AM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.1	T.Flow	: 0.2ml/min
Sample Name :	Bradykinin	B.conc	: 50%H <sub>2</sub> O/50%MeOH
Cat No. :	RP11129		
Lot No. :	PE0631812		
Theoretical MW :	1060.21		
Observed MW :	1060.6		