


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/13/18
Time Processed : 13:50:34
Injection Volume : 0.2
Sample Name : MUC5AC
Cat. No. : RP20290
Lot. No. : PE3721805
Theoretical MW : 1501.62
Observed MW : 1501.60

Interface : ESI
Nebulizing Gas Flow : 1.5 L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH