


# Mass Spectrum



Sample Information  
User : Gary  
Month-Day Processed : 01/10/17  
Time Processed : 1:00:29  
Injection Volume : 1  
Sample Name : Enfuvirtide (T-20)  
Cat No. : RPI7307  
Lot No. : PE3331612  
Theoretical MW : 4491.88  
Observed MW : 4492.20

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH