


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 06/22/18
Time Processed : 10:53:04
Injection Volume : 0.2
Sample Name : HIV-1 TAT Protein Peptide
Cat No. : RP10129
Lot No. : PE0021806
Theoretical MW : 1559.90
Observed MW : 1559.85

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH