


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/23/17
Time Processed : 18:09:28
Injection Volume : 0.2
Sample Name : Brain Natriuretic Peptide(BNP)(1-32),human
Cat .No. : RP11119
Lot. No. : PE4671712
Theoretical MW : 3464.10
Observed MW : 3464.20

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50%H₂O/50%MeOH