


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 11/26/17  
Time Processed : 15:57:04  
Injection Volume : 0.2  
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human  
Cat No. : RP11110  
Lot No. : PE7021710  
Theoretical MW : 2197.60  
Observed MW : 2198.25

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

## Interface Bias

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH