


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 10/20/17
Time Processed : 11:29:21
Injection Volume : 0.2
Sample Name : HIV p17 Gag (77-85)
Cat. No. : RP17316
Lot. No. : PE7051710
Theoretical MW : 981.11
Observed MW : 980.55

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH