

Sample Name :Glucagon-Like Peptide (GLP) II, human
Cat No. :RP10769
Lot No. :E29185211
Time Processed :9:39:03 AM
Month-Day-Year Processed :09/26/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm


Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Chromatogram



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	18.035	18018	2552	0.674
2	18.345	9911	705	0.371
3	18.718	11759	1468	0.440
4	19.080	16202	1214	0.606
5	19.379	8397	850	0.314
6	19.595	14872	1963	0.556
7	19.960	2543585	336230	95.158
8	20.344	50278	7062	1.881
Total		2673022	352043	100.000