


Mass Spectrum



Sample Information

Month-Day Processed :	09/25/21	Interface	:ESI	Equipment	: ZJ21010035
Time Processed :	17:11:39	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Injection Volume :	0.2	CDL Temp	:250	Drying Gas Flow	:5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human	Block Temp	:200	T.Flow	:0.2 ml/min
Cat No. :	RP10769			B.conc	:50%H2O/50%MeOH
Lot No. :	E29185211				
Theoretical MW :	3766.20				
Observed MW :	3765.6				