


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 10/25/18
 Time Processed : 8:25:44
 Injection Volume : 0.2
 Sample Name : TAT peptide
 Cat No. : RP20256
 Lot No. : PE3411810
 Theoretical MW : 1621.91
 Observed MW : 1622.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH