


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	01/18/19	Interface Bias	: +4.5kV
Time Processed :	11:07:32 AM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.2	T.Flow	: 0.2ml/min
Sample Name :	Histatin 5	B.conc	: 50% $\text{H}_2\text{O}$ /50% $\text{MeOH}$
Cat No. :	RP11240		
Lot No. :	PE4281901		
Theoretical MW :	3036.29		
Observed MW :	3036.5		