


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/21/19
Time Processed : 09:00:10 AM
Injection Volume : 0.3
Sample Name : HA Peptide
Cat. No. : RP11735
Lot. No. : E0974819K
Theoretical MW : 1102.15
Observed MW : 1101.8

Equipment : GK11010007
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200