


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/14/16
Time Processed : 15:54:02
Injection Volume : 0.2
Sample Name : Somatostatin
Cat No. : RP10226
Lot No. : PE5061611
Theoretical MW : 1639.90
Observed MW : 1639.90

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH