


Mass Spectrum



Sample Information
Month-Day Processed : 11/02/20
Time Processed : 12:11:49 PM
Injection Volume : 0.3
Sample Name : Magainin 2
Cat No : RP11232
Lot No : E8973120J
Theoretical MW : 2466.90
Observed MW : 2466.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH