


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 02/09/20
Time Processed : 03:38:15 PM
Injection Volume : 0.3
Sample Name : BNP(1-45), mouse
Cat No. : RP20365
Lot No. : E6203020A
Theoretical MW : 4919.60
Observed MW : 4919.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH