


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 06/26/18
Time Processed : 13:27:53
Injection Volume : 0.2
Sample Name : LL37 (Human)
Cat No. : RP13323
Lot No. : PE2761806
Theoretical MW : 4493.28
Observed MW : 4494.40

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Drying Gas Flow : +4.5kV
T.Flow : 5 L/min
B.conc : 0.2ml/min
:50%H2O/50%MeOH