


Mass Spectrum



Sample Information
Month-Day Processed : 06/12/21
Time Processed : 20:33:28
Injection Volume : 0.2
Sample Name : [Glu]-Fibrinopeptide B
Cat. No. : RP20439
Lot. No. : E991121F
Theoretical MW : 1570.60
Observed MW : 1570.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH