


Mass Spectrum



Sample Information

Month-Day Processed :	12/24/21	Equipment	: GK11010007
Time Processed :	12:35:25	Interface Bias	: +3.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	β-Amyloid (25-35)	T.Flow	: 0.2 ml/min
Cat No. :	RP10008	B.conc	: 50%H ₂ O/50%MeOH
Lot No. :	E6790221L		
Theoretical MW :	1060.27		
Observed MW :	1060.6		