


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 09/26/19  
Time Processed : 05:47:00 PM  
Injection Volume : 0.3  
Sample Name : Fibrinopeptide B, human  
Cat No. : RP10836  
Lot No. : E4777819I  
Theoretical MW : 1552.56  
Observed MW : 1551.8

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : -3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH