

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. :RP11110
 Lot No. :E6797021J
 Time Processed : 17:31:01
 Month-Day-Year Processed : 11/12/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: GK11010017

mV



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	12.909	6043	464	0.167
2	13.983	15784	2112	0.437
3	14.136	2467	354	0.068
4	14.479	14446	2731	0.400
5	14.617	3438732	498014	95.138
6	14.977	115562	15070	3.197
7	15.104	13026	2441	0.360
8	15.426	3320	528	0.092
9	15.722	2715	458	0.075
10	15.920	2379	451	0.066
Total		3614474	522623	100.000