


Mass Spectrum



Sample Information
 Month-Day Processed : 11/11/21
 Time Processed : 9:20:07 PM
 Injection Volume : 0.4
 Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. : RP11110
 Lot No. : E6797021J
 Theoretical MW : 2197.60
 Observed MW : 2197.2

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200
 Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH