


Mass Spectrum



Sample Information

Month-Day Processed :	09/18/21	Equipment	: ZJ21010035
Time Processed :	7:59:55	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	Exendin-4	T.Flow	: 0.2 ml/min
Cat No. :	RP10874	B.conc	: 50% H ₂ O/50% MeOH
Lot No. :	E8305621I		
Theoretical MW :	4186.66		
Observed MW :	4186.0		