


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/29/20
Time Processed : 08:07:42 AM
Injection Volume : 0.3
Sample Name : Tetanus toxin (830-843)
Cat No. : RP19980
Lot No. : E3438620E
Theoretical MW : 1611.84
Observed MW : 1611.2

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH