


Version 1.1

Update: 9/28/2019

Certificate of Analysis

Product Information

Product Name:	Human CD8 Antibody (M8A16) [FITC], mAb, Mouse
Cat. No.:	V03005-25T
Quantity:	25 tests
Presentation:	PBS, pH 7.4, containing gelatin and 0.03% ProClin300.
Storage:	The antibody should be stored at 4°C and protected from prolonged exposure to light. Do not freeze.
Purification:	The antibody was purified by affinity chromatography and conjugated with FITC under optimal conditions. The solution is free of unconjugated FITC.
Lot No.:	PF200602
Concentration:	40 µg/mL
Date of Manufacturing:	06/29/2020
Expiry Date:	06/29/2021
Quality control:	Flow cytometric analysis of antibody surface-stained cells.



Human peripheral blood lymphocytes were stained with or without Human CD8 Antibody (M8A16) [FITC].

Date: 06/30/2020

Approved by: *Xianli Yang*

Quality Department

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.