

Version: 01

Update: 01/13/2021

Progesterone Antibody (5F2HC)

Cat. No.	Name	Clone
V09201	Progesterone Antibody (5F2HC)	5F2HC
Specificity	Progesterone	
Isotype	Recombinant human IgG1	
Production	Cultured <i>in vitro</i> under conditions free from animal-derived components	
Purification	Protein A/G affinity column	
Formulation	50 mM Na-citrate, 150 mM NaCl, pH 7.0, containing 0.03% Proclin 300	
Storage	For up to three months at 2-8°C or for up to three years at -25°C~-15°C. Avoid repeated freezing and thawing.	
Application	FIA and CLIA	
Background	The progesterone is a steroid hormone which is secreted by the female reproductive systems. Its concentration is associated with the development and regression of the corpus luteum. The progesterone is used in evaluation of ovulation status and assessment of the luteal phase.	
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.	

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.