

Version: 01  
Update: 11/08/2021

## Testosterone Antibodies

Cat. No.	Name	Clone
V10601	Testosterone Antibody (7B7HC)	7B7HC
V10602	Testosterone Antibody (8A4HC)	8A4HC

<b>Specificity</b>	Testosterone
<b>Isotype</b>	Recombinant human IgG1 for mAbs 7B7HC and 8A4HC
<b>Production</b>	Cultured <i>in vitro</i> under conditions free from animal-derived components
<b>Purification</b>	Protein A/G affinity column
<b>Formulation</b>	50 mM Na-citrate, 150 mM NaCl, pH 7.0, containing 0.03% Proclin 300
<b>Storage</b>	For up to three months at 2-8°C or for up to three years at -25°C~-15°C. Avoid repeated freezing and thawing.
<b>Application</b>	CLIA
<b>Background</b>	Testosterone is an important sex hormone in males. It is synthesized in the body from cholesterol. Testosterone is critical for the development of male reproductive tissues. It also regulates sexual development and red blood cell production.
<b>Note</b>	GenScript can customize this product per customer's request including product size, buffer components, etc.

*For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.*