

Rev04
 Update: Aug,20,2025

DATASHEET


CD3E/CD3 epsilon, His, Cynomolgus

Cat. No.: Z04010


Product Introduction

Species	Cynomolgus
Protein Construction	CD3E/CD3 epsilon (Gln22-Asp117) Accession # Q95LI5 His N-termC-term
Purity	> 90% as determined by Bis-Tris PAGE
Endotoxin Level	Less than 1EU per µg by the LAL method.
Biological Activity	Measured by its binding ability in a functional ELISA. Test result was comparable to standard batch.
Expression System	HEK293
Theoretical Molecular Weight	11.7 kDa
Apparent Molecular Weight	Due to glycosylation, the protein migrates to 14-16 kDa based on Bis-Tris PAGE result.
Formulation	Lyophilized from a 0.22 µm filtered solution in PBS, (pH 7.4).
Reconstitution	Centrifuge the tube before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Storage & Stability	Upon receiving, the product remains stable for 6 months at -20°C or below. Upon reconstitution, the product should be stable for 3 months at -80°C. Avoid repeated freeze-thaw cycles. Recommend to aliquot the protein into smaller quantities for optimal storage.


Examples



CD3E/CD3 epsilon, His, Cynomolgus on Bis-Tris PAGE under reduced condition. The purity is greater than 95%.



Immobilized CD3E/CD3 epsilon, His, Cynomolgus, His Tag at 0.5 μg/ml (100 μL/well) on the plate. Dose response curve for Anti-CD3 Antibody, hFc Tag with the EC50 of 6.2 ng/ml determined by ELISA.



The purity of CD3E/CD3 epsilon, His, Cynomolgus is greater than 95% as determined by SEC-HPLC.

Background

Target Background : CD3E, is a single-pass type I membrane protein. CD3 (cluster of differentiation 3) T cell co-receptor helps to activate both the cytotoxic T cell (CD8 naive T cells) and also T helper cells (CD4 naive T cells). It consists of a protein complex and is composed of four distinct chains. In mammals, the complex contains a CD3γ chain, a CD3δ chain, and two CD3ε chains.

Synonyms : CD3e; CD3E; T3E; FLJ18683; TCRE; CD3-epsilon; IMD18

For research use only. Not intended for human or animal clinical trials, therapeutic or diagnostic use.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China